



Agilent Columns & Supplies keep you on the **fast track** to productivity



Getting started with any new technology can be quite a challenge. Agilent Analyzers and Application Kits help you spend *less time* on set-up and configuration, and *more time* generating the highest quality results from new technologies.

Using decades of experience, customer feedback and collaboration with labs from around the world, Agilent has changed the focus of the installation from meeting general hardware specifications to meeting your specific application needs.



That's why everything included in your Analyzer and Application Kit is designed to help you:

- **Save time and money.** Your application methodology, column, supplies, and checkout standards mixture are included so you won't waste time or money purchasing individual components.
- **Minimize downtime and disruption.** Each Analyzer or Application Kit is configured for your application, includes the hardware and software components needed for optimization and is shipped fully assembled so you get to work fast.
- **Streamline preparation.** Your Analyzer is factory configured and tested, and the checkout report is included at shipping. Simply perform a duplicate on-site checkout and compare results to ensure that your application is ready to go.
- **Get smooth operation right out of the box.** Unlike other providers, Agilent can custom-tailor its Analyzers to your unique needs, and all Analyzers include concise, application-focused video training tools. This means you're ensured excellent results from first use through the entire 10-year life of your instrument.

Agilent would like to continue making your analyses easier – using the same high-quality Agilent columns and supplies that were provided in your Analyzer or Application Kit.

Please see other side for reorder information for GC and GC/MS columns and supplies.

Our measure is your success.



Easy reordering guide: columns and supplies

GC Columns

Application Kit	Description	Part Number
All Pesticides	Agilent J&W HP-5MS Ultra Inert GC Column, 15 m, 0.25 mm, 0.25 µm	19091S-431UI
Japanese Positive List Pesticides	Agilent J&W DB-5MS Ultra Inert GC Column, 30 m, 0.25 mm, 0.25 µm	122-5532UI
Toxicology	Agilent J&W DB-5MS Ultra Inert GC Column, 15 m, 0.25 mm, 0.25 µm	122-5512UI
Oral Fluids	Agilent J&W DB-17MS GC Column, 15 m, 0.25 mm, 0.25 µm	122-4712
Oral Fluids	Agilent J&W DB-1MS GC Column, 15 m, 0.25 mm, 0.25 µm	122-0112
Oral Fluids	Fused Silica Tubing, Deactivated, 0.1 mm, 1 m	160-2635-1

GC Supplies

Application Kit	Description	Part Number
All Pesticides and Toxicology	DB-5MS UI Restrictor 5 m, 0.15 mm, 0.15 µm	165-6626
All Kits	Liner, Helix double taper, Deactivated	5188-5398
All Kits	Gold Plated Inlet Seal with Washer	5188-5367
All Kits	Bleed and Temperature Optimized Septa, 11 mm, 50/pk	5183-4757
All Kits	Certified non-stick liner O-ring, 10/pk	5188-5365
All Kits	Siltite Ferrules, 0.1-0.25 mm Column, 10/pk	5188-5361
All Kits	Blue Line Autosampler Syringe – 10 µL tapered, fixed, Fitted Plunger	G4513-80204
All Kits	Blue Line Autosampler Syringe – 50 µL straight, fixed, PTFE-tipped Plunger	G4513-80221

Standards

Application Kit	Description	Part Number
All Pesticides	Pesticide Locking Standard 3 x 1 mL	5190-1441
All Pesticides	Pesticide Analyzer Check Out Solution	5190-0468
Toxicology	Toxicology Check Out Standard	5190-0471
Oral Fluids	Derivatized THC sample	G1540-85010

Sample Preparation

Application Kit	Description	Part Number
All Pesticides	QuEChERS EN method sample pack	5982-7000
All Pesticides	QuEChERS AOAC method sample pack	5982-7005

For more information regarding the specific products received with your Analyzer/Application Kit please refer to the Quick Start Guide for your system.

Visit www.agilent.com/chem/contactus
to contact your local Agilent Representative
or Agilent Authorized Distributor

